

# Soft Gelatin Capsule Storage Trays - Deeps



800108



624008



805108



Top to bottom  
 887008 lid  
 870008 3" deep  
 875008 4 1/2" deep  
 880008 6" deep  
 870148 dolly



894008

After encapsulation, softgel capsules are soft and sticky due to high moisture content. While the use of a tumble dryer removes a great deal of moisture, additional drying is often needed. MFG Tray Company trays are specifically designed to accommodate the capsules during this period.

- The reinforced composite material provides a smooth surface so capsules aren't marked.
- This material offers excellent dimensional stability and load carrying capacity without creep, distortion, or warp.
- A secure grip provides easy lifting with no sharp edges.

## Deeps - For in-process bulk storage

Model	Description	Outside Top Dimension			Inside Bottom Dimension		
		Length	Width	Height	Length	Width	Depth
624008	Stacking Tray	31	25	2 3/4	28 1/8	23 1/8	2 5/8
800108	Stacking Box	31 3/4	19 1/8	4 1/2	29 3/4	18	4
805108	Stacking Box	30 3/8	15 7/8	2 3/4	29 5/8	15 1/8	2 5/8
870008	Stacking Box	25 3/4	17 3/4	3	23 7/8	15 7/8	2 7/8
875008	Stacking Box	25 3/4	17 3/4	4 1/2	23 7/8	15 7/8	4 3/8
880008	Stacking Box	25 3/4	17 3/4	6	23 7/8	15 7/8	5 7/8
887008	Lid for 870, 875, 880	26 3/4	18 3/4	6			
894008	Trays with drop sides	35 3/4	19 3/4	3 3/4	34	17 7/8	3 1/2