

Soft Gelatin Capsule Drying Trays

Stacking Trays that allow ventilation for drying, curing or aging of softgels, tablets, capsules, powders, teas, and drugs.



Top to bottom:

600208
600008
600108
600308



- Made of thermoset reinforced composite material, MFG Drying Trays offer excellent dimensional stability and load carrying capacity without creep, distortion or warp, in temperatures from -60° F to 250° F. (continuous) or +300° F. (intermittent).
- The sides, or ends, or sides and ends, are dropped so air can circulate over and around the contents to speed drying.
- NSF listed under Standard No. 2, and all ingredients used in the manufacture of our products are listed in US FDA Code of Federal Regulations, Title 21, of the Federal Register, Section 177.2420 acceptable "for repeated used in contact with food".
- These products are used in the drying or cooling process of a number of different products and industries, some of which include candy, pharmaceuticals, minerals, pet food, polymers, and textiles.

630101



633001



641008



642008

Model	Description	Outside Top Dimension			Inside Bottom Dimension		
		Length	Width	Height	Length	Width	Depth
600008	Tray, with drop ends	29 7/8	23 7/8	1 1/2	28 7/8	22 7/8	1 3/8
600108	Tray, with drop sides	29 7/8	23 7/8	2	28 7/8	22 7/8	1 7/8
600208	Tray, with drop sides	23 7/8	14 7/8	1 3/8	22 7/8	13 7/8	1 1/4
600308	Tray, with drop end	29 7/8	23 7/8	2	28 7/8	22 7/8	1 7/8
600408	Tray, with drop end	29 7/8	23 7/8	1	28 7/8	23 7/8	7/8
630101	Tray with drop ends	24 1/4	24 1/4	2 1/4	23 1/4	23 1/4	2 1/8
633001	Tray with drop ends	32 3/4	24 3/4	1 1/2	31 5/8	22 1/2	1 3/8
641008	Tray with drop end	32	24	1 3/4	30	22 1/8	1 5/8
642008	Tray with drop sides & ends	29 7/8	23 7/8	1 1/2	28 7/8	22 7/8	1 3/8